

# KWD-808I

## ELECTRO ACUPUNCTURE STIMULATOR

### Standard accessories :

- + Main machine: 01 unit;
- + Power jack: 01 pcs
- + Fastener
- + Cable stand : 01 pcs
- + White stickers : 2 pair
- + Output wire: 05 group
- + Acupoint detecting pen: 01 piece

### Specifications:

- Detecting Body Acupoint additionally
- KWD-808I impulse electrotherapy can output five kinds of waveform: continuous wave, discontinuous wave, compressional wave, fluctuational wave and start/stop wave.
- The frequency of impulse of the basic wave is from 1.2Hz to 55Hz and the width of it is 0.6ms.
- Single phase main supply 220V / 50Hz or internal 9V battery
- Input power: less than 20 VA.
- Six output channel , stimulate 12 acupoint at the same time
- Adjustable pulse intensity
- Timing function: The timing range is from 0~60 min $\pm$  10%
- Size: 220 x 170 x 75 mm;
- Weight : 2 kg
- Package: 10 boxes / carton



**Wujin Greatwall Medical Device Co., Ltd**

Chengxi Industrial concentrated Area,Hutang Town,Wujin District,Changzhou

Tel:0519-86593745

Fax:0519-86593745