

Elbow CPM



- Progressive mode : Relieves pain by gradual adjustment to the angle of rehabilitating movement therapy.
- Oscillation mode : Maximizes its kinesiatrics effect by repeating exercise at each end- angle.
 - Added wrist rotation mode.
 - Allows to select mode : Uniform / Accelerated mode.
 - Easy operation with TFT LCD(option) or Dedicated LCD.

INFORMATION

Model	ARTUS-501E
Power Supply	110Vac ~ 240Vac, 50/60 Hz

FUNCTIONS

Speed	1 ~ 5	level.
Timer	1 ~ 99	min.
Count	1 ~ 99	
Flexion/Extension Pause	0 ~ 9	sec.

TECHNICAL SPECIFICATION

Flexion / Extension	0° ~ 135°	
Internal / External Rotation	- 90° ~ 50°	
Synchronized abduction and Rotation	Elbow	0° ~ 135°
	Wrist Rotation	- 90° ~ 40°
Speed (Elbow)	60°/min ~ 200°/min	
Speed (Wrist Rotation)	75°/min ~ 200°/min	
Weight	33kg	
Dimensions (WxDxH)	706 × 647 × 835 (mm)	