

IASO^(HBOT)

HYPERBARIC OXYGEN CHAMBER



Mild Hyperbaric Oxygen Therapy · HB-501 · HB-36 · HB-301


SAMBO
VENTEC

Mild Hyperbaric Oxygen Therapy


Henry's law states that the amount of gas dissolved in a liquid or tissue is proportional to the partial pressure of that gas in contact with the liquid or tissue. This is the basis for increased tissue oxygen tensions with HBO treatment.

IASO HB-501


Actual time display & recordable



DISPLAY




SETTING



RECORD





① IASO HB-501

② IASO HB-501 Controller

Specifications

Model	IASO HB-501
Dimension	765mm Ø x 2,200mm
Material	Galvanized steel 2.0mm
Pressure capacity	1.5ATA / 7.5 psi
Safety	De-pressure valve, and relief valve
Control	Micro process
Compressor	120 l/min x300kpa
Cooler	¼ hp

Controlloer

SETTING		Operation time	every minutes
		Pressure	0-50 KPA
		Air cooler	℃
		Ventilation	intermittent
DISPLAY		Air pressure	KPA
		Oxygen concentrate	%
		Temperature	℃
		Operation time	min
		CO2 concentration	ppm
RECORD		Operation time	every minutes
		Pressure	KPA / min
		O2 concentration	% / min